



















Size	Round 	Cutting 
2		
1		
0		
2-0		
3-0		
4-0		
5-0		
6-0		





Synthetic Absorbable Surgical Suture Polyamide, Monofilament

Composition ► Polyamide

Construction ► Monofilament

Sterilization ► Ethylene Oxide

Colour ► Black

Applications ► General closure, Skin and Plastic Surgery.

Features ► Being monofilament, Ilon is remarkably smooth and gives excellent knot security. Possesses good tensile strength. Easy removal with no tissue adherence. Ideal for use in skin closures.