

















Nonabsorbable Surgical Suture, Black Braided Silk

- Composition** ▶ **Black Braided Silk**
- Construction** ▶ **Multifilament**
- Coating** ▶ **Wax**
- Sterilization** ▶ **Ethylene Oxide**
- Colour** ▶ **Black**
- Applications** ▶ **General closure, G. I. Tract, Plastic, Skin Surgeries. Widely used as ligature.**
- Features** ▶ **Coated with wax to allow easy passage through tissue, limit absorption and increase ductility. Knot pull strength and needle attachment strength much more than USP limits.**

Size	Round 	Cutting 
2		
1		
0		
2-0		
3-0		
4-0		
5-0		
6-0	